## **Medical Imaging**

## IST 2022-2023

## Lab 2

Consider the formation using a 1<sup>st</sup> generation scanner of a CT image of the modified Shepp-Logan phantom, which approximately models X-ray attenuation in a human head.

Note: for this lab you will need the <a href="mage">scikit-image</a> processing toolbox.

- 1. Generate the modified Shepp-Logan phantom using shepp\_logan\_phantom. Use the function rescale to get a 256x256 dimension.
- 2. Simulate the sinogram obtained by collecting projections covering [0;180]° in steps of 1° (using radon), and the associated reconstructed image using filtered backprojection (using iradon).
- 3. Repeat the simulation in 2., but now covering [0;90[° and [0;360[°, in steps of 1°. Determine what the minimum angle coverage should be, in theory.
- 4. Repeat the simulation in 2., covering [0;180[°, but now in steps of 5°. Observe the effects on the sinogram and the reconstructed image.
- 5. Repeat the simulation in 2., by adding noise to the projection data, considering a maximum sinogram intensity of 255 (using the random.poisson function of the numpy library; note that values need to be converted to uint8 type (8-bit unsigned integers) prior to applying this function).
- 6. Repeat the simulation in 5., by now considering half of the original tube voltage. Compare with 5.
- 7. Repeat the simulation in 6., by replacing the original Ram-Lak filter by modified filters (available in iradon), and explain the results as a function of their different frequency responses.